
CERUS

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INTERCEPT PROCESSING SET for PLASMA

For use with INTERCEPT Illuminator

Each set is wrapped in a tamper-evident package and includes
 (1) one 15ml 6mM amotosalen hydrochloride solution container
 (Formula : Amotosalen HCl 203mg - Natr. chlorid. 924mg - Aqua ad
 inject. ad 100ml),
 (2) one illumination container,
 (3) one compound adsorption device (CAD),
 (4) three INTERCEPT Plasma storage containers.
 The set is sterilized by a combination of steam and radiation.

Indications and Usage

This set is used with an INTERCEPT Illuminator to inactivate a broad spectrum of viruses, bacteria, parasites and leukocytes in plasma. INTERCEPT Plasma is indicated for support of patients requiring plasma transfusions, according to clinical practice guidelines. INTERCEPT Plasma may be stored for 3 months between -18°C and -25°C and for 12 months below -25°C, from the time of collection, according to applicable procedures and regulations.

Contraindications

Use of INTERCEPT Plasma is contraindicated in patients with a history of allergic response to amotosalen or psoralens.

Precautions

Do not use if : tamper-evident package has been opened; signs of deterioration are visible; fluid path closures are loose or not intact; cannulae are broken or there is no fluid in amotosalen solution container.

Do not store above 25°C. Do not vent. Do not freeze. Protect the pack and tubing from sharp objects.

Unused sets in open aluminium foil may be kept 20 days at room temperature by folding and securing open end of aluminium foil. Units removed from the aluminium foil must be used within 24 hours.

Keep set in light-protective package until time of use. Protect from direct sunlight and strong UVA light source.

Use one set for each plasma unit (385-650ml).

All the following conditions must be met for pathogen inactivation:

- The specific volume ranges of the plasma unit must be maintained.
- Plasma mixed with amotosalen must be exposed to UVA light dose from INTERCEPT Illuminator. No other source of UVA light may be used.
- Plasma must be passed through the CAD by gravity flow process after illumination.
- The entire process from collection to freezing must be completed within a timeframe of 8 hours.

This process is designed to be a closed system. Treatment with INTERCEPT Blood System does not replace applicable standards for processing in open and closed systems. If there is a leak in the set during processing, plasma product must be discarded.

Warnings : Amotosalen in contact with skin may result in photosensitisation in the presence of ultraviolet light. If skin exposure occurs, flush exposed skin copiously with water. Sterile connecting device (SCD) should be used according to manufacturer's directions for use.

Notes to Physicians

While laboratory studies of amotosalen processing with UVA light have shown a reduction in levels of certain viruses, bacteria and parasites; there is no pathogen inactivation process that has been shown to eliminate all pathogens.

Neonatal patients who require plasma transfusion during phototherapy for treatment of hyperbilirubinemia should be treated with phototherapy devices that do not emit light less than 425nm to avoid the theoretical potentiation of an interaction between UVA light and psoralen, resulting in erythema.

A maximum daily intake limit of 1300 µg Amotosalen is recommended.

Instructions for Use

Materials Needed : One (1) INTERCEPT Processing Set for Plasma.

Equipment Needed : INTERCEPT Illuminator, Sterile Connecting Device (SCD), Tube Sealer.

Used INTERCEPT Processing sets should be discarded like any used blood containers, as biohazardous waste.

Process Steps

A-Plasma for Processing with the INTERCEPT Blood System

The plasma for processing with the INTERCEPT Blood System shall have a volume of 385-650ml, including anticoagulant, with red bloodcell content less than 4×10^6 RBC/ml.

B- Amotosalen Addition

1. Remove set from package.
2. Weld tubing from plasma container to amotosalen container tubing using SCD.
3. Disassemble set from organizer and remove rubber band.
4. If two plasma units are desired at the end of the INTERCEPT process, heat seal and remove one storage container.
5. Label set containers using appropriate donation identification system (eg. donation number and product code). Ensure label on illumination container is placed «on flap only». See warning in section C. While labeling storage containers, separate them to ensure they do not adhere to one another.
6. Hang plasma container and break both cannulae on amotosalen container, ensuring that containers/components do not come in contact with floor.
7. Allow plasma and amotosalen solution to flow into illumination container (marked by number «1»).
8. Ensure that plasma drains completely from initial plasma container into illumination container.
9. Express air from plasma within the illumination container into the amotosalen container.
10. When air is removed and plasma has fully drained back into the illumination container, mix illumination container thoroughly by gentle agitation to ensure complete mixing of amotosalen and plasma.
11. Express a small amount of plasma and amotosalen mixture into tubing, filling at least 4cm of tubing. Clamp tubing.
12. Seal tubing between illumination container and amotosalen container within the 4cm.
13. Remove and discard initial plasma container, amotosalen container and excess tubing.

C- Illumination

Illuminate plasma.

Refer to INTERCEPT Illuminator Operator's Manual for complete illumination instructions for use. **Warning :** All plasma, both in illumination container and tubing, must be within large compartment of illuminator tray in order for inactivation to occur. The process requires unimpeded light transmission through tray and illumination container with plasma. No labels or other material should be on this area. Tray must be clean. Labels should be placed on illumination container flap only. Illumination container should lay flat.

D- Processing with Compound Adsorption Device (CAD)

1. Hang illumination container, allowing CAD to hang freely, with storage containers kept in an inverted position.
 2. Close clamp on bypass line; ensure all other clamps are open.
 3. Break cannula on illumination container and allow plasma to flow through CAD (marked with the number «2») into storage containers (marked with the number «3»).
 4. Once plasma has emptied from the illumination container and passed through the CAD, close clamp below CAD.
 5. Hang CAD together with illumination container.
 6. Remove storage containers from tab on CAD and allow them to hang ports up.
 7. Completely express the air from one storage container and close the clamp. Repeat air expression with the second storage container if required. Ensure all air is expressed into one final storage container.
 8. Open the clamp on the CAD bypass line and completely express air from the final storage container into the illumination container. Close the clamp on the bypass line and open the clamps on the storage containers.
 9. Open the clamp below the CAD, allowing plasma to drain into the storage containers.
 10. Ensure that the storage container tubing contains plasma and no air. Close clamps on storage containers.
 11. Re-distribute plasma volume between storage containers, if necessary.
 12. Ensure appropriate donor identification is attached to each filled plasma storage container.
 13. Disconnect storage containers from set by heat-sealing, allowing sufficient tubing length for segments.
 14. Discard CAD and illumination container.
- The INTERCEPT Plasma process is now complete.
15. Seal tubing as appropriate for making segments as needed.
 16. Follow internal procedures for freezing plasma.

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- Cerus, INTERCEPT, et INTERCEPT Blood System sont des marques de Cerus Corporation.