
CERUS

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INTERCEPT PROCESSING SET for LARGE VOLUME PLATELET UNITS

For use with INTERCEPT Illuminator

Each set is wrapped in a tamper-evident package and includes one 17.5ml 3mM amotosalen hydrochloride solution container (Formula : Amotosalen HCl 101mg - Natr. chlorid. 924mg - Aqua ad inject. ad 100ml), one illumination container, one container with Compound Adsorption Device (CAD), one INTERCEPT platelet storage container. The set is sterilized by a combination of steam and radiation.

Indications and Usage

This set is used with an INTERCEPT Illuminator to inactivate a broad spectrum of viruses, bacteria and leukocytes in platelets. INTERCEPT platelets are indicated for support of patients requiring platelet transfusions, according to clinical practice guidelines. INTERCEPT Platelets may be stored up to 7 days from time of collection, at 20-24°C with continuous agitation. INTERCEPT Platelets stored up to 7 days have been shown to adequately prevent and control bleeding. Any extension of platelet storage time from current blood center limits should be evaluated per Directive [2004/33/EC] and validated according to local blood bank procedures.

Contraindications

Use of INTERCEPT platelets is contraindicated in patients with a history of allergic response to amotosalen or psoralens.

Precautions

Do not use if : tamper-evident package has been opened; signs of deterioration are visible; fluid path closures are loose or not intact; cannulae are broken or there is no fluid in amotosalen solution container.

Do not store above 25°C. Do not vent. Do not freeze. Protect the pack and tubing from sharp objects.

Unused sets in open aluminium foil may be kept 20 days at room temperature by folding and securing open end of aluminium foil. Units removed from the aluminium foil must be used within 8 hours.

Keep set in light-protective package until time of use. Protect from direct sunlight and strong UVA light source.

Use one set for each platelet dose.

All the following conditions must be met for pathogen inactivation:

- Platelets should be prepared with InterSol Solution for platelet re-suspension in the volume range according to **Table 1**.
- Specified number of platelets and red blood cells must be maintained according to **Table 1**.
- Platelets mixed with amotosalen must be exposed to UVA light dose from INTERCEPT Illuminator. No other source of UVA light may be used.
- Exposure to UVA light must be completed by end of day after blood collection (end of Day 1).
- Platelets must be agitated in CAD container for a minimum of 6 hours and a maximum of 16 hours after illumination.

Warnings : Amotosalen in contact with skin may result in photosensitisation in the presence of ultraviolet light. If skin exposure occurs, flush exposed skin copiously with water. Sterile connecting device (SCD) should be used according to manufacturer's directions for use.

This process is designed to be a closed system. Treatment with INTERCEPT Blood System does not replace applicable standards for processing in open and closed systems. If there is a leak in the set during processing, platelet product must be discarded.

Notes to Physicians

While laboratory studies of amotosalen processing with UVA light have shown a reduction in levels of certain viruses and bacteria, there is no pathogen inactivation process that has been shown to eliminate all pathogens.

Neonatal patients who require platelet transfusion during phototherapy for treatment of hyperbilirubinemia should be treated with phototherapy devices that do not emit light less than 425nm to avoid the theoretical potentiation of an interaction between UVA light and psoralen, resulting in erythema.

Directions for Use

Materials Needed : One (1) INTERCEPT Processing Set for Large Volume Platelet Units, for each platelet dose (2.5-7.0 x 10¹¹).

Equipment Needed : INTERCEPT Illuminator, Sterile Connecting Device (SCD), Tube Sealer, Flatbed Agitator.

Laboratory temperature should be controlled to meet applicable regulations for platelet processing.

Used INTERCEPT sets should be discarded like any used blood containers, as biohazardous waste.

Process Steps

A- Preparation of Platelets

Platelets should be prepared with InterSol solution. Platelets are suspended in approximately 35% plasma and 65% InterSol solution. Platelet products within the following ranges have been shown to be acceptable for processing with INTERCEPT Blood System.

Table 1

| Platelet Count | Final Volume | Plasma content | RBC Content |
|------------------------------|--------------|----------------|-----------------------------|
| 2.5 – 7.0 x 10 ¹¹ | 300-420ml | 32-47% | <4 x 10 ⁶ RBC/ml |

B- Amotosalen Addition to Platelets

1. Remove set from package. Unwrap only illumination container from organizer.
2. Weld tubing from platelet pooling container to amotosalen container tubing using SCD.
3. Label set containers using appropriate donation identification system.
4. Hang platelets and break both cannulae on amotosalen container.
5. Allow platelets and amotosalen solution to flow into illumination container, marked by number "1".
6. Ensure that platelets are completely transferred to illumination container by expressing air from platelet and amotosalen container into illumination container.
7. When platelets are in illumination container, mix thoroughly by gentle agitation to ensure complete mixing of amotosalen and platelets.
8. Express air from platelets into amotosalen container.
9. Express a small amount of platelet and amotosalen mixture into tubing, filling about 4cm of tubing.
10. Seal tubing between illumination container and amotosalen container, so that tubing is no longer than approximately 4cm from illumination container inlet port. **Warning :** Tubing must be held within large compartment of illumination tray.
11. Remove and discard empty platelet and amotosalen containers.

C- Illumination

Illuminate platelets.

Refer to INTERCEPT Illuminator Operator's Manual for complete illumination directions for use. **Warning :** Platelets in illumination container must be within large compartment of illuminator tray in order for inactivation to occur. The process requires unimpeded light transmission through tray and illumination container with platelets. No labels or other material should be on this area. Tray must be clean. Labels should be placed on illumination container flap only. Illumination container should not be folded.

D- Agitation with CAD

Warning : Do not fold or bend CAD.

1. Unwrap set from organizer.
2. Hang platelets, break cannula and allow platelets to flow into CAD container, marked by number "2".
3. Express air from CAD container into illumination container.
4. Seal tubing close to inlet port of CAD container.
5. Remove and discard empty illumination container.
6. Place CAD container on flatbed agitator for at least 6 hours, but no more than 16 hours.

E- Transfer to Storage Container

1. Remove platelets from agitator and hang platelets.
2. Close clamp on platelet sampling pouch.
3. Break cannula and allow platelets to flow into storage container, marked by number "3".
4. Express air from storage container into CAD container.
5. Seal tubing close to inlet port of storage container.
6. Remove and discard empty CAD container. Place platelets on platelet agitator.

The INTERCEPT Platelet process is now complete.

F- Sampling Platelet Product (optional)

1. Mix INTERCEPT platelets by gently agitating platelet storage container.
2. Open clamp to platelet sampling pouch and squeeze several times.
3. Allow sampling pouch to fill with platelets. Seal tubing.
4. Remove sample pouch.
5. Transfer sample to appropriate laboratory tube immediately.

Cerus, INTERCEPT, and INTERCEPT Blood System are trademarks of Cerus Corporation

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Cerus, INTERCEPT, et INTERCEPT Blood System sont des marques de Cerus Corporation.