



May, 2008

Re: **INTERCEPT Blood System™
Plasma Processing Set (Product Codes INT310X)
Maximum Starting Plasma Volume Expansion**

Dear Valued Customer,

CERUS is committed to providing you with quality products and services. This letter is to inform you of a change that is being implemented in our INTERCEPT Blood System product line. We are pleased to advise you of an improvement to the approved labeling for the INTERCEPT Plasma processing sets which has been made in response to your feedback.

Details of Improvement

The maximum plasma volume for processing with the INTERCEPT Blood System has increased from a 635ml to 650ml. As a result, the new range for plasma processing with the INTERCEPT system is 385ml - 650ml, including anticoagulant (excluding the 15ml Amotosalen added during processing).

This change has been through conformity assessment by the Notified Body, TÜV SÜD Product Service. The updated INTERCEPT Processing Set for Plasma Instructions For Use reflecting this improvement will be available in Q2 2008.

Existing inventory of disposable sets packaged with the current Instructions For Use (IFU) will continue to be distributed until such inventory is depleted; however, if you'd like, you may contact Cerus BV Customer Services to request a copy of the new Instructions For Use. Upon availability, and after such request, the new IFU will be sent to you.

We trust you will find this improvement useful and are available for any questions or inquiries.

Contact information

If you have any questions concerning this change, please contact CERUS BV Customer Services:

Stationsstraat 79-D
3811 MH Amersfoort
The Netherlands

Email: customer_services@cerus.com
Phone: +31 33 496 0600
Fax: +31 33 496 0606

With Kind Regards,

Cerus Europe BV